

BinderTutorial_QC

March 6, 2019

1 Binder Tutorial QC Analysis

1.0.1 To begin: Click the top cell and press 'Run' on the toolbar (or shift-enter). Alternatively click Kernel, Restart and Run All.

1.1 Table of Contents:

1. Section ??
2. Section ?? 2.1. Section ?? 2.2. Section ?? 2.3. Section ?? 2.4. Section ??

1. Import Data

1. Import the pandas python module (<https://pandas.pydata.org/>).
2. Import the excel sheet "Data" from excel file "data.xlsx" into a data frame called "data".
3. Display the numbe of rows and column.
4. Display the fist 10 rows at the top (head) of the data frame.

```
In [10]: import pandas as pd
```

```
data = pd.read_excel('data.xlsx', sheet_name='Data')

print("Data Table: {} rows & {} columns".format(*data.shape))
display(data.head(10))

peak = pd.read_excel('data.xlsx', sheet_name='Peak')
print("Peak Table: {} rows & {} columns".format(*peak.shape))
display(peak.head(10))
```

Data Table: 91 rows & 3087 columns

	Order	SampleType	QC	M1	M2	M3 \
0	1	QC	1	1.837062e+08	1.296528e+08	3.961182e+07
1	2	QC	1	2.030571e+08	1.178121e+08	6.361418e+07
2	3	Sample	0	1.285728e+08	2.726444e+08	4.611692e+07
3	4	Sample	0	1.491128e+08	1.955086e+08	5.212802e+07
4	5	Sample	0	1.590268e+08	3.780061e+08	4.238071e+07
5	6	Sample	0	1.875314e+08	7.153474e+07	5.007826e+07
6	7	QC	1	1.983585e+08	1.204062e+08	5.190127e+07

7	8	Sample	0	1.460268e+08	2.274148e+08	5.872350e+07
8	9	Sample	0	1.381372e+08	1.269111e+08	5.276344e+07
9	10	Sample	0	1.601298e+08	4.347327e+08	4.883041e+07

	M4	M5	M6	M7	...	\
0	3.732578e+07	5.339865e+06	1.031822e+08	9.409682e+06	...	
1	6.682747e+07	5.572493e+06	9.366510e+07	9.941482e+06	...	
2	4.717621e+07	3.178726e+06	8.077425e+07	6.924520e+06	...	
3	5.114760e+07	3.428193e+06	8.211311e+07	7.416907e+06	...	
4	4.282051e+07	3.186309e+06	1.000761e+08	7.030797e+06	...	
5	5.036435e+07	4.849348e+06	9.998757e+07	1.039036e+07	...	
6	5.180519e+07	5.348706e+06	8.966852e+07	9.605938e+06	...	
7	5.581725e+07	4.342786e+06	1.059740e+08	1.401639e+07	...	
8	4.936756e+07	3.005090e+06	8.101362e+07	1.484466e+07	...	
9	4.518047e+07	3.326398e+06	7.714394e+07	1.239396e+07	...	

	M3075	M3076	M3077	M3078	M3079	\
0	35153.206002	15735.886504	24928.165047	30134.444665	3035.247672	
1	36701.059662	15590.522963	24776.756320	30250.013996	3052.288747	
2	43120.310719	16378.981747	28627.451582	31297.776297	1738.358685	
3	44663.545552	17429.459487	25703.703594	26399.809078	846.020446	
4	41770.351723	18467.519040	27371.393974	28374.785933	1959.397898	
5	41616.726360	11426.528395	23313.373795	28582.727888	872.115128	
6	32839.675140	16262.163327	21298.643088	28065.010551	3277.276918	
7	35339.968664	25426.050392	51819.646796	21762.932664	1822.815140	
8	33766.256647	15963.341703	50651.507533	18503.876558	843.925809	
9	31442.161072	22600.725041	51901.700174	16804.403434	1522.101762	

	M3080	M3081	M3082	M3083	M3084
0	18271.101975	20981.783849	3523.418763	23752.440055	17161.698785
1	22116.296519	24543.382677	3389.677457	24262.317947	18949.290754
2	17695.679626	16703.470005	1113.489872	29717.457693	14683.792853
3	15306.546820	15198.095655	1137.752119	26610.868749	11913.223789
4	14904.243032	15625.089253	1101.738536	26496.665356	12750.259715
5	15301.063816	14321.556052	1030.238432	32587.650554	18655.811727
6	20071.616958	18198.098904	3047.831055	24997.154651	14494.628663
7	23457.954245	23651.858710	1210.741328	22834.426860	15771.899526
8	16572.962662	16625.898352	1186.166100	25623.108950	13811.758227
9	19948.743747	22690.058851	1229.664436	22637.770945	15009.816887

[10 rows x 3087 columns]

Peak Table: 3084 rows & 6 columns

Idx	Name	Mol_Weight	RT_minutes	RSD	D_Ratio
0	1	M1	113.05902	1.276	5.012937
					0.28


1	2	M2	203.11564	0.950	9.557013	0.19
2	3	M3	161.10514	1.428	16.085023	0.26
3	4	M4	129.07893	1.400	18.035797	0.49
4	5	M5	161.10510	1.272	19.515359	0.40
5	6	M6	225.94404	5.923	6.510319	0.33
6	7	M7	194.08026	1.169	2.890900	0.03
7	8	M8	115.06345	7.067	4.733385	0.16
8	9	M9	131.09461	6.220	6.138996	0.33
9	10	M10	117.07900	6.971	4.783892	0.32

2. Visualisation

Note: Each cell in the Visualisation Section can be run in any order (provided data is imported in Section 1). ### 2.1. Histogram of RSD

```
In [7]: import matplotlib.pyplot as plt
        %matplotlib inline

fig = plt.figure(figsize=(6,5))
plt.hist(peak.RSD, 50, density=True, facecolor='g', alpha=0.5)
plt.xlabel('RSD', fontsize=15)
plt.show()
```




2.2. Jointplot of RSD vs. D-Ratio

```
In [11]: import seaborn as sns
```

```
sns.jointplot(x=peak.RSD, y=peak.D_Ratio, kind='kde', color="skyblue")
```

```
Out[11]: <seaborn.axisgrid.JointGrid at 0x12a5c2eb8>
```



2.3. PCA score plot of QC vs. Sample

```
In [8]: import numpy as np
        from sklearn.decomposition import PCA
```

```

from sklearn.preprocessing import StandardScaler

import matplotlib.pyplot as plt
from matplotlib.colors import ListedColormap

names = peak['Name']
x = data[names].values
x = np.log(x)
x = StandardScaler().fit_transform(x)

pca = PCA(n_components=2)
scores = pca.fit_transform(x)
label = data['SampleType']

Sample_scores = scores[label == 'Sample',:]
QC_scores = scores[label == 'QC',:]

fig = plt.figure(figsize=(8,8))


h1 = plt.scatter(Sample_scores[:,0],Sample_scores[:,1],edgecolors='Black', facecolors=
h2 = plt.scatter(QC_scores[:,0],QC_scores[:,1], edgecolors='Black', facecolors='Red',s

plt.legend((h1,h2),('Sample','QC'),fontsize=15)

plt.xlabel('PC1', fontsize=15)
plt.ylabel('PC2', fontsize=15)
plt.title('Quality Control PCA plot',fontsize=20)


plt.show()

```



2.4. Scatter plot of Molecular Weights vs. RT Time (sized by RSD)

```
In [9]: fig = plt.figure(figsize=(20,16))
plt.scatter(peak.Mol_Weight,peak.RT_minutes, s = peak.RSD**2/2,alpha=0.2, edgecolors =
plt.xlabel('Molecular Weight', fontsize=15)
plt.ylabel('RT minutes', fontsize=15)
plt.title('Metabolites Detected (sized by RSD)',fontsize=20)
plt.show()
```



In []: